

DIMAR

PRODUCT LINE

Intensive Care Unit
Hospital Departments

DimAir Line
The Helmet
REVOLUTION



ZIP HELMETS

Single patient Mechanical Ventilation (NPPB) helmet, also suitable for PSV-ASB, BIPAP and CPAP. Equipped with a hermetic zip, they allow quick and easy access to the patient. Convenient and practical.

TECHNICAL FEATURES:

- Pneumatic zip opening*
- Made of biocompatible plastic.
- 2 fastening eyelets on the top of the helmet for an alternative/ additional fixing.
- Symmetrical R/L connections for convenient and easy connection with equipment on either side.
- Sealed access from 3 to 19 mm.
- Soft and flexible armpit straps to prevent the decubitus of the armpits, complete with protection in antibacterial material CoolMax®.
- Connections ISO 22M, 22F.
- Latex and PHT free.



NIMV 6R ZIP

Non-invasive ventilation helmet (NPPB), usable in PSV-ASB, Bilevel and CPAP. 6 ISO connectors placed on the rigid ring for greater Effectiveness, Practicality and Comfort. The connection to the hard ring reduces the response delays of triggers <0.1".



ADVANTAGES:

- Device not in contact with the patient's face.
- Better seal.
- Quick access to the patient's face, also during emergencies.
- Possibility of breastfeeding, feeding and weaning without removing the helmet.
- 6 symmetrical connectors (3 on the right and 3 on the left) to allow positioning of the flow generator/ventilator on both sides of the bed and offer more flexibility in patient mobilization.

CODE	SIZE	RDM	CND/EMDN
600/5586	XXS	2683511	R0301010402
600/5234	XS		
600/5235	S		
600/5236	M		
600/5237	L		
600/5238	XL		
600/5239	XXL		
600/5582	3XL		



COMFORT CPAP ZIP

Helmet for CPAP. Without rigid ring, prevents decubitus to the patient's nape and shoulders. The high comfort allows you to use it for long periods in the ward.

CODE	SIZE	RDM	CND/EMDN
500/5579	XXS	2683503	R0301010401
500/5441	XS		
500/5442	S		
500/5443	M		
500/5444	L		
500/5445	XL		
500/5446	XXL		
500/5568	3XL		



COMFORT NEO-PED ZIP

CPAP helmets for neonate patients from 3 to 5 Kg, and pediatric from 5 to 15 Kg. Without a rigid ring, prevents decubitus at the patient's nape. The retaining pad is made of antibacterial CoolMax® hygienic.

CODE	SIZE	RDM	CND/EMDN
500/5447	NEO	2683503	R0301010401
500/5448	PED		

- In single pack
- Multiple pack of 5
- Possibility to create customized models, versions and packaging on specific request



HMC Premedical S.p.A

Via Bosco 1/3 – 41037 Mirandola – MO – I
Tel: +39 0535 22704 Fax: +39 0535 609546

