

**DIMAR**

PRODUCT LINE

**DimAir Line**

### PEEP VALVES

PEEP VALVES adjustable from 0 to 20 cmH<sub>2</sub>O, suitable for maintaining positive pressure in Respiratory systems.

Single pack. Multiple pack 10 pieces.

Pre-calibrated from 2.5 to 20 cmH<sub>2</sub>O, code 700/6339.



CODE	DESCRIPTION	RDM	CND/EMDN
700/5056	PEEP VALVE 0 - 20 cmH <sub>2</sub> O	2683521	R028001
700/6336	PEEP VALVE 0 - 20 cmH <sub>2</sub> O with spring		
700/5049	PEEP VALVE 0 - 20 cmH <sub>2</sub> O WITH ANTIMICROBIC FILTER		
700/5050	PEEP VALVE 0 - 20 cmH <sub>2</sub> O WITH SPRING AND ANTIMICROBIC FILTER		



### CPAP BALLOON

High compliance device ideal for stabilizing the pressure.

Capacity of 25 liters. Connector 22Flex. Single pack. Multiple pack 10 pieces.

CODE	DESCRIPTION	RDM	CND/EMDN
700/5055	CPAP BALLOON	2683516	R03020102

### MANOMETER

To measure the pressure value in the CPAP system.

Graduated scale from 0 to 20 cmH<sub>2</sub>O.

Single pack. Multiple pack of 10 pieces.



CODE	DESCRIPTION	RDM	CND/EMDN
700/6355	22F HELMET MANOMETER	2683488	R03010180
700/6353	22M HELMET MANOMETER		
700/3342	MASK MANOMETER		

**CPAP ACCESSORIES**



## TWO-WAY VALVE :

**CODE 700/5067**

Anti-choking/asphyxia valve. Opens and closes automatically at the limit pressure value of 2.5 cmH<sub>2</sub>O. The automatic opening of the valve allows fresh air to enter and the gas to exit. Single pack. Multiple pack of 10 pieces.



## ONE-WAY VALVE:

**CODE 700/5058**

One-way safety valve. Prevents risk of asphyxia in case of pressure loss in positive pressure respiratory systems. Single pack. Multiple pack 10 pieces.



## OVER PRESSURE VALVE:

**CODE 700/5057**

Safety valve pre-calibrated to 20 cmH<sub>2</sub>O. Applied on CPAP and NIMV helmets, it prevents the increase of internal pressures to levels above 20 cmH<sub>2</sub>O. Single pack. Multiple pack 10 pieces.



## CVC 4 WAYS EXTENSION:

**CODE 700/3358**

For connection with pneumatic seal to central venous catheter non-invasive ventilation with the helmet. Pneumatic seal. Single pack. Multiple pack 10 pieces.



## AEROSOL KIT:

**CODE 700/5065**

Kit for administration of nebulized drugs. Composed of a mouthpiece, connecting tube, nebulizer and oxygen supply tube. Single pack. Multiple pack of 10 pieces.



**HMC Premedical S.p.A**

Via Bosco 1/3 – 41037 Mirandola – MO – I  
Tel: +39 0535 22704 Fax: +39 0535 609546



0425