

DIMAR

PRODUCT LINE

Intensive Care Unit &
Hospital Departments

100%
MADE IN ITALY

DimAir Line

INTERCHANGEABLE
HEADGEAR

FIFTH
FIXATION POINT

VERSIONS
- SINGLE PATIENT
- RESTERILIZABLE

MORE SPACE FOR
THE NOSE

NASOGASTRIC
TUBE
DEDICATED
PASSAGE

INTERCHANGEABLE
CONNECTORS

CONNECTION
GAS INLET

TOTAL FACE
DIMAXZERO

"DIMAX ZERO" TOTAL FACE MASK

New generation Total Face facial mask for non-invasive mechanical ventilation in PSV, BIPAP and CPAP techniques. The innovative configuration reduces CO₂ rebreathing to zero. Featuring 5 attachment points, it optimizes patient/ventilator synchronization and pressure distribution to the face, reducing leakages and improving comfort. Wide space for prominent noses. The high volume dedicated to the nose and the possibility of replacing the connectors to transform the same mask into Vented or Non Vented, make DIMAX ZERO the ideal interface for non-invasive ventilation.

8 plus

- + Zero rebreathing
- + Connection
- + Comfort
- + Sealing
- + Practicality
- + Flexibility
- + Safety
- + Versions

CO₂ rebreathing
Fluid-dynamic
inside
Total Face
Standard
Mask compared
to DIMAX ZERO



TECHNICAL FEATURES:

- Zero dead space.
- 5 fixation points.
- More space for the nose.
- Dedicated access for NasoGastric tubes.
- Single patient and resterelizable.
- Interchangeable Connectors:
 - 22F Non Vented
 - 22M Vented



ADVANTAGES:

- Zero CO₂ rebreathing.
- Better patient/ventilator synchronism.
- Improved pneumatic seal.
- More comfort for prominent noses.
- Lower headgear traction.
- Reduced risk of contact injury.
- Possibility of passage for nasogastric tube and fibrobronchoscopy, without leakage and interruption of ventilation.
- Single Patient and resterelizable in autoclave at 134°C - gas plasma – ETO versions.

REFERENCES:

SINGLE PATIENT**	SIZE	RDM	CND/EMDN
700/6919	XS/PED	2683312	R0301010202
700/6828	SMALL		
700/6829	MEDIUM		
700/6830	LARGE		
700/6831	EXTRA LARGE		
700/7226	EXTRA EXTRA LARGE		

** In single pack

*** In single pack

RESTERELIZABLE***	SIZE	RDM	CND/EMDN
700/6920	XS/PED	2683333	R0301010202
700/6824	SMALL		
700/6825	MEDIUM		
700/6826	LARGE		
700/6827	EXTRA LARGE		
700/7225	EXTRA EXTRA LARGE		

Multipack of 5



The 5th fastening point allows greater grip on the patient's forehead



Dedicated access for NasoGastric tubes



Connection "NonVented" for bi-tube ventilators



Configuration for CPAP with venturi generator "Easy Vent"



Configuration with bi-tube ventilators



Configuration with single tube ventilators or CPAP generator



DIMAX ZERO XS/Pediatric

The images and contents reproduced are indicative and non-binding

5018/01 FIRST ISSUE FOR NEW CODING